

# Cytek Biosciences DxP Athena<sup>TM</sup> Flow Cytometry System

Built by Cytek from the ground up, the DxP Athena<sup>™</sup> flow cytometry system incorporates Cytek's proven DxP technology enabling the capability to resolve dim populations in a multicolor format. Choose from configurations with up to 3 lasers and 13 fluorescent detectors. Each configuration offers consistent powerful performance to meet your application requirements at an affordable price.

DxP Athena benefits include:

- Intuitive software
- A small footprint
- High sensitivity and resolution for resolving dim populations
- Built-in maintenance capability
- Upgradeability
- Optional Cytek 96-Well Automated Micro-Sampler (AMS) and FlowKart™ Fluid Management System
- World-class Cytek Service

For more information, visit Cytek at www.cytekbio.com.



## **Optics**

#### **Excitation Optics**

Optical platform

Allows up to 3 lasers. Fixed optical assembly with three spatially separated laser beams.

Lasers

407nm: 50mW 488nm: 45mW 637nm: 80mW

Beam geometry

Prismatic expander and achromatic

spherical lens.

Optical Efficiency
Power loss at flow cell: <20% of specified

laser power

**Emission Optics** 

Optical coupling

Fused silica cuvette coupled to emission lens by refractive index matching optical gel for optimum collection efficiency.

Background rejection

Obscuration blade and slit for minimizing unwanted laser radiation at the detector.

Forward scatter detector and filter High-performance solid-state silicon detector with 488nm bandpass filter.

Side scatter detector

Photomultiplier and 488nm bandpass filter.

**Photomultiplier Tubes** 

High efficiency Tri-alkali Photomultipliers

Fluorescence detectors and filters DxP Athena 13-color Configuration shown. Other configurations available.

BluFL1: 530/30 (FITC, Alexa Fluor® 488)

BluFL2: 575/30 (PE)

BluFL3: 615/25 (PE-Texas Red®) BluFL4: 695/40 (PerCP-Cy5.5)

BluFL5: 780/60 (PE-Cy7) RedFL1: 661/16 (APC)

RedFL2: 710/50 (Alexa Flour® 700)

RedFL3: 780/60 (APC-Cy7)

VioFL1: 450/50 (Brilliant Violet™ 421) VioFL2: 525/50 (Brilliant Violet™ 510) VioFL3: 615/25 (Brilliant Violet™ 605) VioFL4: 710/50 (Brilliant Violet™ 711) VioFL5: 780/60 (Brilliant Violet™ 785)

### **Fluidics**

#### Sample Flow Rates

Front panel keypad provides four modes: Run, Standby, Prime, and Clean.

Three preset flow rates:

LO: 12 µL/min
MED: 35 µL/min

HI: 60 µL/min

Maintenance

Onboard CLEAN mode for monthly

maintenance.

Standard Fluidic Reservoirs One 4-L sheath container and one 4-L waste container provided.

# DxP Technology Performance

Fluorescence Sensitivity
Molecules of Equivalent Fluorescence (MEFL),
using Q&b Method\*

| Fluor        | Min. Q | Max. b | **Typical R<br>Value |
|--------------|--------|--------|----------------------|
| FITC         | 0.007  | 1200   | 382                  |
| PE           | 0.05   | 1000   | 300                  |
| PerCP -Cy5.5 | 0.003  | 600    | 336                  |
| PE-CY7       | 0.001  | 2000   | 758                  |
| APC          | 0.015  | 500    | 161                  |
| APC-CY7      | 0.005  | 7500   | 1161                 |
| BV421        | 0.01   | 3750   | 818                  |
| BV510        | 0.01   | 1500   | 1261                 |

<sup>\*</sup>Q measures optical efficiency, b measures background, and R (resolution limit) measures the number of dye molecules required to resolve a dim population from noise

Fluorescence Sensitivity Threshold FITC: 50 molecules of equivalent soluble flurochrome (MEFL-FITC)

PE: 30 molecules of equivalent soluble flurochrome (MEFL-PE)

\*FITC and PE measurements performed using SPHERO Rainbow Calibration Particle (RCP-30-5A)

Fluorescence Linearity Delivers doublet/singlet ratio of 1.95–2.05 for CEN stained with PI and excited with the 488nm blue laser.

Forward and Side Scatter Sensitivity Enables separation of fixed platelets from noise

Forward and Side Scatter Resolution Performance is optimized for resolving lymphocytes, monocytes, and granulocytes.

Side Scatter Resolution Capable to resolve 0.5µm beads from noise

Fluorescence Resolution 18 bit 5 log decades

Data Acquisition Rate 7,500 events/s with beads.

<sup>\*\*</sup>Average R value across 4 systems. MEFL required to be 2 standard deviations above noise. Assumes no compensation applied.



# Data Management

FlowJo<sup>™</sup> Collector's Edition 7.5.110 or later

Our acquisition interface can be adapted to any application.

Real time spillover matrix for viewing live compensated data.

Acquisition templates include hardware, gate, spillover, layout and statistical settings.

Save workspaces and use during multiple acquisition sessions.

Cytek AMS software Version 1.0.4

Fast and Easy Setup Experiment plate mapping, including stain names can be setup on stand alone computer and saved as a template for future use.

Block Inspector Change the FlowJo Collector's Edition acquisition settings on a per block basis.

Well ID Stamping
Well ID in the file name confirms data file to
Well ID relationship.

## Workstation

Operating System
Windows® 7 32-Bit Professional

Processor Intel Quad Core processor, 3.0 GHz

4GB (1 x 4GB), 16000 MHz DDR3

Hard Drive 500GB SATA 3.0Gb/s

DVD Drive 16x DVD+/- RW, SATA

Video Processor HD Graphics GMA 4600

Monitor 24" LCD

## **Options**

96-Well Auto-sampling

Cytek's automatic micro-sampler (AMS) can be mounted directly under the sample injection tube reducing dead volume and increasing throughput.

Auto-mode

Full 96-well plate in 40 minutes (15 sec. acquisition/well and 10 sec. wash/well) less than 1% carryover.

High throughput mode Full 96-well plate in 15 minutes. (9 sec. acquisition/well) less than 3% carryover.





# Installation Requirements

*Dimensions (W x D X H)* Sensor module (without AMS) 55 x 52.4 x 57.8 cm (21.7 x 20.6 x 22.8 in)

Weight 45.4 kg (100lbs)

Computer 48 x 41 x 58 cm (19 x 16 x 23 in)

Recommended workspace (W x D x H) 180 x 91 x 132 cm (71 x 36 x 52 in)

#### Room Requirements

Power 100-240V, 50/60 Hz, 2A max

Heat dissipation 450 watts with all solid-state lasers

Temperature 16–29°C (60–85°F)

Humidity 10% to 90% relative non-condensing

Air filtering
No excessive dust and smoke

Lighting
Optics and detectors shielded from room

# Regulatory Status

For Research Use Only. Not for use in diagnostic or therapeutic procedures except in countries where the regulatory approval has been obtained from the local regulatory authorities. Please check with your local sales representatives before placing your orders.



Cytek Biosciences, Inc. 46107 Landing Parkway Fremont, CA 94538 (877) 92-CYTEK sales@cytekbio.com